

A black and white close-up portrait of a smiling woman with long, dark dreadlocks. She is looking directly at the camera with a bright, joyful expression, showing her teeth. The background is dark and out of focus.

In-Kone[®]

A fresh take on
«**Bone Level**»
implantology

The logo for Global D, featuring a stylized orange circular icon with a white, curved, 'G'-like shape inside.

Global D

Contents.

01.

In-Kone[®],
a fresh take on the "Bone Level" philosophy _____ p7

02.

Technical characteristics promoting
triple integration _____ p10

03.

A comprehensive range of implants
to suit all clinical situations _____ p28


04.

ULTIMATE
Surgical protocol _____ p36



Our philosophy: **Perio-implantology** in essence

Global D is a French manufacturer with more than 30 years' experience in the design of dental implants, with a constant concern for the **preservation of the biological space**. This approach enables us to propose a unique and proven range of therapeutic solutions, **developed to promote the bone, mechanical and periodontal integration of implant-based prosthetic restorations**.



TRIPLE INTEGRATION, a key **SUCCESS** factor for your implant-based restorations

Mechanical integration

Immobility and integrity of the prosthetic assembly during masticatory movements.^{(2) (3)}

Periodontal integration

Barrier against the risks of leakage and contamination at the implant-prosthesis interface.^{(1) (4)}

Osseous integration

Osseointegration of the implant and preservation of the peri-implant bone plates.⁽⁵⁾

These three types of integration are closely related and must be combined to **optimise the overall result** of the restoration, both aesthetically and functionally.


At Global D, we are particularly attentive **to this overall vision, which guarantees a successful and durable treatment plan**; it serves as our reference system for the development of our products.

Triple integration is enhanced by the right balance between a sufficiently large implantable tissue space and an aesthetic and functional prosthetic plan.

A **global** approach focused
on **aesthetic & functional** durability

Analysis and preparation of the implantable space

1. The volume and quality of the hard and soft tissues must be able to ensure the **mechanical anchorage**, immune **response** and **aesthetic appearance** of the implant-based restoration. An adequately prepared tissue environment is a prerequisite for the dental implant to perform its function biomechanically in a reliable manner.




Respect of positioning

2. The three-dimensional positioning of the implant is determined by both the **tissue environment**, which defines the "implantable space", and the rehabilitation plan, which defines the "**prosthetic space**". The implant emergence axis must fulfil the requirements specific to each of these two spaces to ensure the stability of the restoration in the long term.

Analysis and design of the prosthetic space

3. The prosthetic restoration is designed to ensure optimal preservation of **the health and stability of the peri-implant tissue capital**. Quality of assembly seal and stability, passivity of frameworks, choice of prosthetic materials, occlusal adjustment and clinical monitoring of restorations are all key factors for successful prosthetic integration in the long term.



01.

In-Kone®,
*a fresh take on the
"Bone Level" philosophy*

Our core concept: one of the first sealed prosthesis-implant connections and a deliberate subcrestal position enable the bone plates to be widened, thus providing natural and long-lasting support to the mucosa.

“

A DELIBERATE SUBCRESTAL POSITION

In-Kone® is a two-piece dental implant that can be used in one or two surgical steps, as well as for aesthetic or immediate loading procedures. It falls into the so-called “Bone Level” category of implants, with the particularity of enabling subcrestal positioning.

It has a unique roughened chamfered collar and internal connection via a friction cone, two characteristics that enable the implant shoulder to be positioned **from 1.5 to 2mm below the cortical bone**.

MEASURABLE CLINICAL BENEFITS*

Widening of the **bone plates**

The subcrestal position recommended for the In-Kone® implant promotes **the rapid formation of a cervical blood clot** and **releases tension in the cortical bone** to allow the latter to grow horizontally onto the implant collar.

* Doctors Patrick MOHENG, Philippe ROCHE-POGGI, Romain CASTRO, Pierre MOHENG, Laboratoire Guillaume PÉNARANDA Implant Global D In-Kone® SA² - Survival rate at 5 years - Implant magazine 2019

In-Kone®

Survival rate
At 1 year: **99.4%**
At 5 years: **99.10%**
Longitudinal study of 1,788 implants*

Mucosa volume

The enlarged bone plate around the In-Kone® implant provides **natural support for the peri-implant soft tissue**. The “tulip-shaped” profile of the associated prosthetic components then promotes **crimping of the restorations** with a **thick mucosal seal**.⁽⁴⁾



02.


Technical characteristics
promoting **triple integration**



PERIODONTAL INTEGRATION

A **stable** implant-prosthesis connection

With the aim of preserving the tissues, the In-Kone® prosthesis connection (**8° friction cone type**) is designed to propose the most **stable** connection possible, both **mechanically** and in terms of its **seal**. The connection stability prevents the risks of micro-movements caused by the occlusal forces and metal deformation, two phenomena that can cause bacterial leakage at the interface⁽¹⁾. Finally, this stability means that it is possible to propose prosthetic components with a concave design at the connective tier of the peri-implant mucosa⁽¹⁾.



Horizontal section view



Vertical section view

Respect **of the anatomy**
of each patient




A WIDE RANGE OF HEALING SCREWS


suitable for all types of gums and clinical cases to enable modulation of the prosthetic emergence required.

As well as different transgingival heights (h) and diameters, the range of healing screws proposes flat or high head formats, to enable optimal management of all clinical situations⁽⁴⁾.

The transmucosal profile of the healing screws is perfectly **homothetic** to that of the final components, but with a very slightly larger diameter (+0.4mm) to avoid any tension on the connective fibres when the final prosthesis is put into place (for both high and flat versions), enabling optimal management of **all situations**⁽⁴⁾.



h: transgingival height in mm



TULIP-SHAPED PROSTHETIC COMPONENTS

The range of In-Kone® prosthetic components offers a **broad choice of references** for progressive management of different biological heights. The “tulip-shaped” transmucosal profile of the prosthetic components is designed to ensure an **atraumatic** fit with the prosthetic seating formed by the corresponding healing screw.

The concave emergence prosthetic components of the In-Kone® range are designed to **free the thick connective tissues and encourage their reorganisation**. This promotes vascularisation⁽⁴⁾ and limits the risk of telescoping with the bone edges.


A colour code indicating the selected periodontal height enables simple, immediate matching of the healing screws and prosthetic abutments.




CUSTOMISABLE SOLUTIONS

Global D strives to adapt to the changes in your practices, and therefore proposes a **wide range of prosthetic solutions that can be customised in the laboratory:**

- Direct implant titanium bases for customised zirconia abutments or screw-retained teeth
- Blanks (or pre-milled) to make personalised titanium abutments
- Titanium bases for conical abutments for implant bars and bridges
- Customised frameworks in **Global D Compliance** certified laboratories for validated connections.



OSSEOUS INTEGRATION



An implant collar to promote **bone bonding**

The anatomical collar of the In-Kone® implant, roughened and chamfered at 45°, frees the cortical bone from all mechanical stress during the healing phase. **It promotes bonding with the bone** and crimping of the implant⁽⁵⁾, two factors that facilitate tissue management and the formation of papillae.



Two **complementary** implant shapes

The **UNIVERSAL** profile has a gradual double thread to achieve **more stable primary anchorage** in low density bones. This profile is ideal for **immediate loading** rehabilitation.

2.0 mm pitch 0.8 mm pitch


360° rotation
=
Apical progression **2.0 mm**

360° rotation
=
Apical progression **0.8 mm**

A single **drilling protocol**

The **PRIMO** profile has a regular single thread suitable for surgeons seeking a certain **flexibility** in high density bones.

The two In-Kone® implant shapes, PRIMO & UNIVERSAL, use exactly the **same ULTIMATE drilling protocol**, for site preparation, allowing the surgeon to choose the profile best suited to the bone density, respecting the receiving site and encouraging osseointegration of the implant.




A carefully selected titanium alloy

Our In-Kone® implants are manufactured from **medical grade TA6V ELI titanium alloy**. This alloy enables better **machining precision** than with “pure” titanium. This is necessary to achieve a **precise, sealed** friction cone type connection, whose **mechanical resistance** is sufficient to prevent crushing of the implant-abutment assemblies over time⁽²⁾.

Recognised surface treatment


The SA² surface finish of In-Kone® implants is obtained by sandblasting followed by double etching. The resulting roughness (Ra) is between 1.5 and 2 microns. This is described as moderate roughness by scientific literature and it encourages both **cellular differentiation** and the proliferation of osteoblasts on the surface of the implant⁽³⁾. A study of In-Kone® implants on beagles showed **contact osteogenesis at 3 weeks** and **bones in the remodelling phase** in the implant thread spaces at 12 weeks⁽³⁾.



1.5 to 2 µm



3 weeks



12 weeks

Certified clean commitment

The SA² surface finish has been certified by the German **Clean Implant Foundation**, which guarantees the cleanliness of implants based on random market sampling.

The members of the scientific board of the **Clean Implant Foundation** are internationally-reputed dental implantology researchers⁽³⁾.

www.cleanimplant.com ►



MECHANICAL INTEGRATION




Mechanical stability of the **8° friction cone**

This friction cone offers an **effective solution** for **prosthetic stability** thanks to the homogeneous distribution of stress on the implant-abutment contact surfaces and the **reduction of micro-movements** under the effects of occlusal loads.



Respect of the **connection**

To guarantee the **performance of the abutment assembly** and implants, Global D uses an extremely precise **machining process**. This quality of assembly has repercussions on the periodontal, bone and mechanical aspects, which means that manufacturing tolerances must not be altered during the restoration phases. Machining or reworking of direct implant connections are not authorised.



To preserve the quality of the prosthetic connection, Global D supplies prosthesis laboratories with parts such as **titanium bases or blanks with ready-machined connections**. The customisable part is then handled by the **prosthesis laboratories**.

With the aim of sharing skills for the dentistry profession, Global D also invites the laboratories that can meet the equipment requirements to become **Global D Compliance** partners.

This scheme enables them to obtain programmes and tools to make implant bridges and bars on straight or angled conical abutments, with no titanium base, according to a recognised quality charter.


Planned **prosthesis-guided** surgery



The **three-dimensional positioning** of implants is a key factor of the durability of implant-related restorations, so Global D implants are found in most of the planning software products available on the market.


Prosthesis-guided implant planning, associated with piloted or guided surgery, has become an integral part of Global D's digital therapeutic offer.

Global D proposes a full **piloted surgery kit** and a full guided surgery kit, based on application of the ULTIMATE protocol.



03.

A comprehensive range **of implants**
to suit all clinical situations



In-Kone® |W|


In-Kone® |N|

In-Kone® |S|

In-Kone® |W|

In-Kone® |N|

CORE CONCEPT




The original so-called standard (ST) prosthetic platform, the **very backbone** of the In-Kone® philosophy, is proposed with 4 implant diameters to suit the vast majority of your clinical cases.

This platform is essential for your routine actions, and is ideal for **full arch restoration** and the **treatment of anterior sections (single or multiple cases) and posterior sections (multiple cases)**.

In-Kone® |ST|

4 implant diameters,
1 standard connection



ST platform In-Kone® implants are available in **UNIVERSAL** and **PRIMO** versions.



After operation

Ideal for
aesthetic work

The In-Kone® implant was originally designed to treat most situations encountered in dental implantology.



Before operation




Before operation

The external profile of the **In-Kone® UNIVERSAL** implant is a cylindrical cone, equipped with a gradual double thread.

This shape is particularly suitable for **post-extraction implantation with immediate loading**.

The possibility of gradual management of the emergence profile, thanks to its tulip-shaped prosthetic components, makes this implant ideal for interventions with an aesthetic purpose.



After operation

The narrow alternative, opting for *finesse*


The In-Kone® concept is also available with an implant diameter of 3 mm **for restoration of small mesiodistal spaces of the maxillary and mandibular lateral incisors.**

In this context, the implant is equipped with a 5° friction cone and a smaller prosthetic platform, identified as Narrow (NR). **The 3.0 implant (NR)** therefore has its own range of prosthetic components.

In the case of multiple replacement of mandibular incisors, make sure each tooth is replaced with an implant and avoid the use of pontics.




Ø 2.3 mm



Ø 3 mm

In-Kone® |NR|

|WD|
In-Kone®



The wide alternative, opting for *resilience*

In response to the evolution of clinical and prosthetic practices (manufacture of machined one-piece crowns, immediate treatment of single posterior edentulous spaces, etc.) and the changing patient population (younger patients, development of bruxism, etc.), Global D has extended its In-Kone® range with a larger prosthetic platform, indicated Wide (WD).

This option is a solution for clinical cases requiring up to 60% higher mechanical resistance of the implant-abutment assembly, without compromising the original philosophy of the In-Kone® implant.

WD platform In-Kone® implants are particularly recommended for the **rehabilitation of single and multiple posterior cases.**



04.

The **ULTIMATE** surgical protocol

A single surgical kit for
all In-Kone® implants

Our ULTIMATE surgical protocol stands out with its homothetic drills and its compatibility with all our implants, regardless of their shape (UNIVERSAL and PRIMO) or platform (ST, NR, WD).


Latest generation drills

The drills have a double cutting edge and a helical groove to ensure **the evacuation of bone chips**.



A simple, user-friendly protocol

Obvious markers for easy reading to **enable management of the different drilling depths** as the drills move back and forth during irrigation. A colour code for each diameter marked on the drills and pins in the kit enables **immediate reading** of the drilling sequences to be implemented.



Modularity of the ***final drilling diameter***

The very gradual incrementation of different drill diameters allows the final drilling diameter to be adjusted to **different bone densities**.

Each reamer preserves the homothetic shape of the preparation, regardless of the final drilling diameter.

FOR MORE INFORMATION ON THE FRESH TAKE ON THE BONE LEVEL PHILOSOPHY



www.globald.com



FIND OUT ABOUT OUR EXCLUSIVE OVERALL APPROACH TO THE OROFACIAL AREA



BIBLIOGRAPHY

(1) Prosthetic gap and microbial leakage

- Kano et al. A classification system to measure the implant-abutment microgap. Int. J. Oral Maxillof. Impl. 2007; 22:879-885
- Dibart et al. In vitro evaluation of the implant-abutment bacterial seal: The lockin taper system. Int. J. Oral Maxillof. Impl. 2005;20:732-737
- Gross M et al. Microleakage at the abutment-implant interface of osseointegrated implants. A comparative study. Int. J. Oral Maxillof. Impl. 1999; 14:94-100
- Quirynen M et al. Microbial penetration along the implant components of the branemark system. An in vitro study. Clin Oral Impl. Res. 1994; 5:239-244
- Yuya Sasada. Implant-Abutment Connections: A Review of Biologic Consequences and Peri-implantitis Implications. Int. J. Oral Maxillof. Impl. 1304 Volume 32, Number 6, 2017
- Hermann et al. Influence of the size of the microgap on crestal bone changes around titanium implants. A histometric evaluation of unloaded non submerged dental implant. A radiographic in the canine mandible. J periodontol. 2001; 72: 1372-1383.
- King GN et al. Influence of the size of the microgap on crestal bone level in non-submerged dental implant. A radiographic study in the canine mandible. J periodontol. 2002; 73: 1111-1117
- Hermann et al. Factors influencing the preservation of the peri-implant marginal bone. Implant Dent. 2007; 16: 165-175.
- Ericsson I et al. Different types of inflammatory reactions in peri-implants soft-tissues. J Clin Periodontol 1995; 22: 255-261
- Ericsson I et al. Radiographical and histological characteristics of submerged and non-submerged titanium implant. An experimental study in the labrador dog. Clin Oral Impl. Res. 1996; 7: 20-26
- Hermann et al. Biologic width around one and two-piece titanium implants. Clin. Oral Impl. Res. 2001;12: 559-571

- Broggini N. et al. Persistent acute inflammation at the implant-abutment interface. J. Dent. Res. 2003; 82: 232-237
- Broggini N. et al. Peri implant inflammation defined by the implant-abutment interface. J. Dent. Res. 2006; 85: 473-478
- Mishra SK Microleakage at the Different Implant Abutment Interface: A Systematic Review, J. Clin. Diagn. Res. 2017 Jun; 11(6): ZE10-ZE15

(2) Conical interface stress

- Luca Fiorillo Finite Element Method and Von Mises Investigation on Bone Response to Dynamic Stress with a Novel Conical Dental Implant Connection. Hindawi BioMed Res. Intern. Volume 2020, Article ID 2976067, 13 pages.
- Tanja Rack et al. An In Vitro pilot study of abutment stability during loading in new and fatigue loaded conical dental implant using synchrotron-based radiography. Int. J. Oral Maxillof. Impl. VOL 28, Number 1 2013.
- Zipprich H et al. A new experimental design for bacterial microleakage investigation at implant-abutment interface: an in vitro study. Int. J. Oral Maxillof. Impl. 2016; 31:37-44
- Zipprich H et al. The micromechanical behavior of implant-abutment connections under a dynamic load protocol. Clin. Impl. Dent Rel. Res; 2018 Oct;20(5):814-823

(3) Implant surface finish

- Al Nicolas-Silvente. Influence of the Titanium Implant Surface Treatment on the Surface Roughness and Chemical Composition. Mat. Impl. Dent. And Reg. Med. DOI: 10.3390/ma13020314
- Albrektsson T, Wennerberg A Oral implant surfaces: Part 1--review focusing on topographic and chemical properties of different surfaces and in vivo responses to them. Int J Prosthodont. 2004 Sep-Oct; 17(5):536-43.

- Elias CN, Meirelles L. Review Improving osseointegration of dental implants. Expert Rev Med Devices. 2010 Mar; 7(2):241-56.
- Wennerberg A, Albrektsson T. Effects of titanium surface topography on bone integration: a systematic review. Clin Oral Implants Res. 2009 Sep; 20 Suppl 4():172-84.
- Wennerberg A, Albrektsson T. On implant surfaces: a review of current knowledge and opinions. Int J Oral Maxillofac Implants. 2010 Jan-Feb; 25(1):63-74.
- Dirk U. Duddeck. On the Cleanliness of Different Oral Implant Systems: A Pilot Study. J. Clin. Med. 2019, 8, 1280; doi:10.3390/jcm8091280

(4) Peri-implant mucosa

- Nozawa T. et al. Biologic height-width ratio of the buccal supra-implant mucosa. Eur. J. Esthet. Dent. 2006, 1, 208-214.
- Farronato D. Ratio between Height and Thickness of the Buccal Tissues: A Pilot Study on 32 Single Implants. Dent. J. 2019, 7, 40; doi:10.3390/dj7020040
- G. Finelle. Peri-implant soft tissue and marginal bone adaptation on implant with non-matching healing abutments: micro-CT analysis

(5) Bone behaviour

- Alberto Monje. The Critical Peri-implant Buccal Bone Wall Thickness Revisited: An Experimental Study in the Beagle Dog The Int. J. Oral Maxillof. Impl. 2019 DOI: 10.11607/jomi.7657
- D. Boskaya. Evaluation of load transfer characteristics of five different implants in compact bone at different load levels by finite element analysis. J. Prosthet. Dent 2004 Dec; 92(6):523-30
- Dr C. Bolle. Soft-tissue and marginal bone-adaptation on platform switched implants with a Morse cone connection: a histomorphometric study in dogs. Inter. J. Period. Rest. Dent. 2016, 36 (2), pp.221-228.

OUR PUBLICATIONS



Early Periimplant Tissue Healing on 1 Piece Implants with a concave transmucosal design: A histomorphometric study in dogs
C. Bolle & coll., Implant dentistry / Vol. 24, Number 5, 2015

Soft tissues and marginal bone adaptation on platform switched implant with a Morse cone connexion: a histomorphometric study in dogs
C. Bolle & coll., The International of Periodontics & Restorative Dentistry. Vol 36 Number 2, 2016

One Piece Implants with Smooth Concave Neck to Enhance Soft Tissue Development and Preserve Marginal Bone Levels: A Retrospective Study with 1 to 6 Year Follow Up.
JP Axiotis & coll., Hindawi Bio Med Research International Volume 2018, Article ID 2908484, 7 pages

Posterior jaws rehabilitated with partial prosthesis supported by 4 mm x 4 mm or by longer implants: four month post loading results from a multicenter randomized controlled trial.
M. Esposito & coll., Euro J oral implantol. 2015; 8 (3): 110

Posterior jaws rehabilitated with partial prosthesis supported by 4 mm x 4 mm or by longer implants: One year post loading results from a multicenter randomized controlled trial.
P. Felice & coll., Euro J oral implantol. 2016; 9 (1): 35-45

Standard drilling versus ultra sonic implant site preparation: A clinical study at 4 weeks after insertion of conical implants.
C. Makary & coll., Implant Dentistry, Vol. 26 Number 4, 2017

Four year post loading results of full arch rehabilitation with immediate placement and immediate loading implants: a retrospective controlled study.

A. Simonpieri & coll., Quintessence International Volume 48, Number 4, April 2017
Finite Element Method and Von Mises Investigation on Bone Response to Dynamic Stress with a Novel Conical Dental Implant Connection.

Luca Fiorillo & coll., Hindawi BioMed Research International Volume 2020, Article ID 2976067, 13 pages

These products are Class I, IIa and IIb Medical Devices and therefore bear CE markings pursuant to Directive 93/42/CEE. These devices are not subject to French social security reimbursement. The medical devices shown may not be available for sale in all countries. For any additional information, please contact the Global D sales department. Before use. If you have any doubts, please contact the Global D sales department. Certain instructions may be available in electronic format. In this case, a QR code and URL will be shown on the device label. However, they are also available upon request and at no additional cost within 7 days. Send all requests to: quality@globald.com

Find all our solutions
on www.globald.com

00IC-CONCEPT_IK-EN-Ed1-22



Products **manufactured in FRANCE**

Global D
ZI de Sacuny
118 Avenue Marcel Mérieux
Brignais - France
Tel. +33 (0)4 78 56 97 00
Fax +33 (0)4 78 56 01 63
www.globald.com

FIELDS OF APPLICATION



Implantology
Oral surgery
Pre-implant surgery
Orthognathic surgery
Reconstructive surgery
Facial trauma surgery
Cancer surgery
Craniomaxillofacial surgery
Orthodontics
Training